

Startup Procedure

• Log into FOM

• Before lighting the plasma, check...

1. Argon pressure (in hallway) ~90psi
2. Water circulator (in hallway) ~20°C
3. Exhaust (inside instrument)
4. Sample introduction- Clamp down pump tubing
 Empty waste
 Fill uptake container with **fresh** di-H₂O
 Fill Autosampler rinse bottle with 3% HNO₃

• Starting ICP-MS:


1. Turn on Autosampler (on lab cart) **Check rinse container: 2-3% acid**
2. Open MS- Thermo Plasmalab Service (on taskbar near clock)
3. Start PlasmaLab (on desktop)

• Go to 'Instrument' tab on left toolbar



• Look for 'Vacuum ready'

• Click the 'On' button (change to operate state) & wait until you see 'Operate'



• Check task bar readings:

	Definition	Value	If	Then
Neb	nebulizer back pressure	~2bar	low	leak or broken neb
			high	clog
Fwd	Forward power- minor pg	1412W		Set based on how good rf generator is
Ref	Reflected power from plasma	< 5 W	high	automatically adjust load & tune
Load	rf power	not zero or 255		
Tune	rf power	not zero or 255		
Exp	Expansion pressure	2mbar	>2.5	hole in the cone is too big or sample cone is in wrong Or vac pump
Ana	Analyte pressure	2.9x10 ⁻⁷		
Speed	Turbo pump speed	1000Hz		
Load	What pump uses	0.8-1.2A	>3	vac leak

- Let the instrument settle for 10mins before tuning
- (Turn on CCT gas during these 10mins if you are using that mode -it's the 1st CCT-He slider bar in 'Add. Gases')
- (Turn off the CCT gas before tuning)

Autosampler Controls

- Make sure Autosampler is checked in **Configurations** (Cetac 260)

Available Accessories and Devices:


	Use	Model	Device Type	Description
	<input checked="" type="checkbox"/>	Cetac 260	Autosampler	serial port 9, baud rate 9600
	<input type="checkbox"/>	Cetac 500	Autosampler	serial port 1, baud rate 9600
	<input type="checkbox"/>	MicroProbe II	Laser	UP213 Local
	<input type="checkbox"/>	Advantech PCI-1620	Trigger	Using un-filtered inputs
	<input type="checkbox"/>	Cetac 500	Autosampler	serial port 1, baud rate 19200
	<input checked="" type="checkbox"/>	Peri Pump	Peri-pump	operate at 34%

- Choose Accessories window button  (under **Tune** or **Configurations**)
- Choose 'Autosampler' tab
- Initialise the accessory unit 

Go to Rinse station Go to Sample Start/Stop Rinse station pump

Choose Rack Type

Choose Sample position



Peristaltic pump setup Autosampler

Connection: COM 9; 9600 baud Rinse height: 140

Go to rinse at the end of an experiment
 Go to rinse after connecting
 Use Additional rinse

Configured racks:
Rack 1 type 60
Rack 2 type 60

21 (3 * 7 * 160)

Tube Position:
Rack: 1
Row: 1
Column: 1
Height: 144

Additional Rinse:
Rack: 0
Row: 0
Column: 0
Height: 0

OK Cancel Apply

- Click OK when finished (do not un-initialise )

Peri-pump

- Normal pump speed is 25 rpm

Tuning- All Users

Standard Mode

- Wait about 10mins while the ICP warms up (*fresh* di-H₂O uptake)
- Go to **Configuration** tab-
 - Check X Series Default in ‘Configurations’
 - Check Autosampler: Cetac 260
and Peri-Pump in ‘Accessories’
 - Pick latest ‘Instrument Settings’ (top of the list)
- If the instrument was not in ‘X Series Default’ mode (ie. left in CCT)
 - Go to **Tune** tab and click Reset button 
 - **Again, this is ONLY IF it was not in ‘X Series Default’**

1a. Autotuning Lens: (everyday)

- Use Standard Tuning solution (Multi-Element **Tune A**)
 - Click Start the Autotune Wizard button  (under **Tune** tab)
 - Choose:
 - Sequence ‘**XSII Xt- without Co**’ (check desktop)
 - Introduce sample ‘manually’
 - It takes approximately 10 mins
- ★ • ‘View report’ to see if Autotune passed & click ‘Finish’ ★
- If it fails, find me.

Replace uptake tubing into di-H₂O while setting up experiments

- Now is a good time to rinse out the autosampler lines. 
- Refer to the information in the previous section
- Go to Rinse station, and watch your element of interest

Experienced users only:

Standard Mode- Performance & Calibration

2b. Performance Report: (everyday, if you manually checked performance)

- Use Standard Tuning solution (Tune A)
- Click Start the Performance Report Wizard button  (under Tune)

Choose: 'XSII Xt'

Introduce sample 'manually'

It takes approximately 5 mins

- ★•'View report' & scroll all the way to the bottom to see if Performance Report passed & click 'Finish' ★

- If it Fails because: sensitivity or interference problem, Autotune
Mass Calibration fails, Find me

3. Detector Cross Calibration: (do once a month)

- Use Multi-Element Tune F solution

Solution contains a high concentration of low mass elements and low concentration of high mass elements to produce the same number of counts on the detector.

- Click Launch the Instrument Calibration Wizard button  (under Tune)

Choose: 'Detector cross calibration'

Introduce sample 'manually'

It takes approximately 5 mins

- Click 'Finish'

Tune F stock soln \approx 50ppm Dilute 1/1000 TuneF working soln \approx 50ppb
--

4. Mass Calibrate the Quadripole: (do once a month)

- Use Standard Tuning solution (Tune A) takes about 3mins

- Click arrow next to Instrument Calibration button  (under Tune)

- Choose Mass Calibrate the Quadripole

Choose Introduce sample 'manually'

- 'View report' & click 'Finish'

- Go into file 'Calibrations>Mass Calibration' – Highlight Standard

Click off outliers:

Peak Width- 0.75 +/- 0.1AMU (0.65-0.85AMU)

Error- Above 0.1 click off (right axis)

Right click on graph to look at one or the other

Experienced users only:

Standard Mode Tune

1b. Manually Check Performance: (everyday, if you are experienced)

- Wait about 15mins while the ICP warms up (di-H₂O uptake)
- Go to **Configuration** tab-
 - Check X Series Default in 'Configurations'
 - Check Autosampler and Peri-Pump in 'Accessories'
 - Pick latest 'Instrument Settings' (top of the list)
- Go to **Tune** tab and click Reset button  (if above was not checked)
- Put uptake line into Standard 10ppb Tuning solution (Multi-Element **Tune A**) & restart Real Time Display 
- Look for Interference's (CeO & Ba++) Ratio < 3%
Stability < 10%

Sensitivity:

Analyte	Instrument Specification	Standard
7Li	>60,000	140,000
115In	>400,000	1,000,000
238U	>800,000	1,700,000
CeO/Ce	<0.02	0.015
Ba++/Ba+	<0.03	0.02

Co spec
150,000

- **If:** Sensit. is low/ Interfer. are high
Stability is high (>10)
CeO/Ce is high
Ba++ is high
Li sensitivity is low
U sensitivity is low
- **Then:** run autotune
change uptake tubing
lower Nebulizer manually
raise Nebulizer manually
raise Focus manually
lower Focus manually
- Save the settings after manually adjusting lenses 
- Run the Performance Report

Lens voltages:

Extraction	-100 to -300
Lens 1	-800 to -1400
Lens 2	-50 to -130
Focus	5 to 15
D1	-20 to -60
D2	-120 to -180
Pole Bias	0 to 0.5
Hexapole Bias	-4 to -3
Neb	0.84 mL/min
Lens 3	-180 to -200
DA	-30

Replace uptake tubing
into di-H₂O while setting
up experiments

CCT Users Only:

CCT Mode Tune

1. Run CCT1 Gas for 5-10 mins during instrument warm-up

2. Tune 1st in Standard Mode

3a. Tune CCT Mode (option 1):

- Use Standard Tuning solution (**Tune A**) without As
- Go to **Configurations** tab and check 'CCT 1 Default'

Don't click Re-use to keep Std tunes settings ()

- Click Autotune button  (under **Tune**)

Choose:

'Quick CCT Tune'

Introduce sample 'manually'

It takes approximately 5 mins

- 'View report.' If Report passed, click 'Finish'
- If it doesn't pass, usually you can purge the CCT gas longer & it will pass. Find me if it doesn't.

CCT specs:

Analyte	Instrument Spec
78Se	< 20 cps
115In	> 100,000 cps

Lens voltages:

Extraction	-100 to -300
Lens 1	-800 to -1400
Lens 2	-50 to -130
Focus	0 to -5
D1	-20 to -60
D2	-120 to -180
* Pole Bias	-14 (-9)
* Hexapole Bias	-17 (-6)
Neb	0.84 mL/min
Lens 3	-180 to -200
* DA	-20
* CCT -He/H	2.8 to 4.0

*different from Std Mode

Replace uptake tubing into di-
H₂O while setting up
experiments

Setting up an Experiment

When to use Std Mode:

- Use this setting for most analytes
- Always use for elements below mass 23 (CCT is not good)
- Watch out for:
 - Analytes with IE (ionization energy) greater than 9 may be a problem. You may need a longer dwell time for better results. Default is 10ms so try 50ms.
 - Analytes with low second IE, like Column I & II elements, may produce doubly charged ions (use cool plasma)
 - In a high sodium solution, NaAr interferes with Cu-63. Use Cu-65.

When to use CCT Mode:


- To analyze elements that have polyatomic interferences.
- Common elements: V, Cr, Fe, As, Se
 - 51V: In a solution of HCl, ClO interferes.
 - 52Cr: Watch out for ArC
 - 53Cr: In a solution of HCl, ClO interferes.
 - 56Fe: In a high calcium solution, CaO and ArO interferes.
 - 57Fe: In a high calcium solution, CaOH interferes.
 - 75As: In solution of HCl, ArCl interferes.
 - 77/80/82Se: Watch out for ArCl, Ar₂, ArHArH
- Other possible elements: Mg, Si, P, S, K, Ca, Ni, Co, Cu
 - 24Mg: High carbon soln, watch out for C₂
 - 28Si: Watch out for N₂
 - 31P: Watch out for NOH
 - 34S: Watch for O₂, detection limit is still 60ppb at best.
 - 39/40K: Watch out for ArH, Ar
 - 40Ca: Watch out for Ar
 - 58/60/62Ni: High calcium soln, watch for CaO, CaOH
 - 59Co: High calcium soln, watch out for CaO, CaOH
 - 63Cu: High sodium soln, watch out for ArNa
 - 66Zn: Watch out for ArMg

Setting up an Experiment, cont.

•Check pump speed: should be 25 rpm (Inst. tab, Accessories window )

Starting an Experiment from an Existing Experiment

- Go to ‘Experiment’ tab on the left toolbar
- Create an experiment from an existing experiment
- Locate the file that you want to use as a template



Timings Pop-up

- When it asks ‘Do you want to change the timings?’ Say ‘No’
- This is the Uptake/Washout delay (25 sec under maximum delay)

Skip directly to the Sample List: (check the other tabs for accuracy)

Sample List:

- Rack: 0- Standard Rack (All row 1)
 - 1- Farthest Sample Rack
 - 2- Closest Sample Rack
- Row/Column: Each rack is labeled, columns go right to left
- Show advanced tab: Remove Int. Stds. here if they look wrong in results.
Make a new Int. Std list. Choose this new list for each sample that has a problem in the list.

★ Fully Quantitative Concentrations:

- Fill in the your standard concentrations for each analyte



Connect the Uptake line to the Autosampler:

Go to the rinse station for a couple mins,
if you have not already done so.




To start your run: press the Queue button,
save your file name, and press Append.

Setting up an Experiment, cont.

•Check pump speed: should be 25 rpm (Inst. tab, Accessories window )

•Starting a new ICP experiment:



- Go to 'Experiment' tab on the left toolbar
- Create a new blank experiment (usually choose 'Continuous')
 - Continuous (ICP only)
 - TRA- time resolved analysis (profile for LA & HPLC)
- Select a database (usually choose 'default'):
 - default (standard ICP settings)
 - default cct (Collision Cell ICP settings)
 - multiple mode database 2005 (both std & cct settings)

•Setup:


- Configuration Editor: should already be highlighted correctly & unchecked
 - X Series Default (std mode) or CCT 1 Default (CCT mode)
 - Autosampler & Peri-pump in Instrument tab
 - Latest Instrument settings are highlighted
- Timings: default is 25 sec delay- under 'maximum delay'
 - *Au, Ag, and Pt need 60 sec washout. You should have HCl too. (Check Uptake/Washout if you want to monitor an element. Make sure to include a minimum delay if you check this.)
- Analyte:
 - Double click both your analytes and internal standards
 - Single click on an analyte to choose multiple isotopes
 - Elements that have known interferences will automatically turn them on. You will see a blue dot for those elements.
 - Pb: Percentages vary so edit equation to add 1x each isotope
- Acquisition Parameters:
 - Survey: 10 sweeps, 0.6ms dwell, 10 channels, std resolution
 - Main-Peak jump: 100 sweeps, 10ms dwell, 1 channel, std res.
 - To increase resolution on main runs, you could:

Lower sweeps to 50 & raise dwell to 50ms

-For high Na samples, use 5ms dwell and high resolution.

-To reduce interferences in your analyte peak, you could:

Change resolution to 'High' to change peak width



Setting up an Experiment, cont.

•Setup, cont:

•Internal Standards: One by one, choose your Internal Standards from the left list and move to the right list 

-Good internal standards are made to produce 100,000 to 300,000 counts. You may need to increase the concentration of the lower mass elements in the internal standard set to achieve this. Typically 2- **5ppb** are used.

-Technique: Interpolate- applies mass curve using all the internal standards and corrects analyte based on curve. If one int. std starts to drop too low and therefore drift (Ex- low mass elements at the wrong concentration), you can take that std out per sample in the Sample List.

^{115}In - corrects analyte based on In internal std only
This std should be within 60 mass units of analyte.
CCT mode: Int std should have similar IE.

•Isotope Ratio: not usually necessary

•Instrument Calibrations:

•Current Mass Calibration & Detector Cross Calibration are automatically used so you do not need to do anything on this page.

•Calibration Method:

•Method:★None- Internal Standards

Semi-quant- if you want to analyze an element and don't have a standard for it (maybe use in survey run)

Accurate within 20%

Addition- Standard addition

★Fully-quant- Normal calibration method for your analytes.

-When semi-quant is checked on the right, this analyte will be used to generate the semi-quant curve.

-Uncheck Semi-quant for Pb if you use an equation to calculate results (see setup>analyte).

Setting up an Experiment, cont.

•QC Setup:

- This section is only used in environmental labs.

•Sample List:

•Sample List:

- Type: Blank- di-H₂O + acid + Internal Standards
Fully Quant Standard- Multiple ppb standards + acid + Internal Standards
Unknown- Your samples + acid + I.S.
- Rack: 0- Standard Rack (All row 1)
1- Farthest Sample Rack
2- Closest Sample Rack
- Row/Column: Each rack is labeled, columns go right to left
- Show advanced tab: Remove Int. Stds. here if they look wrong in results.
Make a new Int. Std list. Choose this new list for each sample that has a problem in the list.

★•Fully Quantitative Concentrations:

- Fill in the your standard concentrations for each analyte

•Connect the Uptake line to the Autosampler:

Go to the rinse station for a couple mins,
if you have not already done so.



To start your run: press the Queue button,
save your file name, and press Append.

Monitoring your Run

Results tab:

- Calibration Data: Check your calibration curves for contaminations
Click on a point in the graph to remove a standard from the curve
Intercept Conc- Detection Limit
- Numerical Results:
 - **Analyte Dilution Conc**- Corrected for internal standards and dilution (specified in sample list)
Watch Internal Standards through out the run ($\approx 0.7-2.0$)
If Internal Standard Correction goes above 2.0, the cones are probably getting clogged & results may be off.
 - **Mass Uncorrected ICPS**- Watch Internal Standard Raw counts (115In $\approx 300,000$ for 5ppb)
Record Raw Indium counts at the beginning and end of run in the Log Book
- Spectra:
 - Check to make sure impurity peaks are not overlapping with your analyte peaks (see Setup>Acq. Parameters>Resolution)
 - When you see a peak half way between mass units, you probably have doubly charged ions from an element at 2x that mass number.
 - The dotted line represents the analog signal and the solid line is the digital. They should match if the detector is calibrated correctly.

Reports

Two options:

- Generate a Report using the PlasmaLab software
Go to Report tab
Choose: Fully Quant Calibrations
Dilution Corrected Concs, etc.
Save or Print
- Generate an Excel file:
Right click on your Numerical Results
Export as csv or copy

Shutdown

•Rinse:

- Rinse in Autosampler rinse station
- Go to Instrument>Tune>Accessories window button 
Go to Peri-pump tab>Initialise 
- Click on Rabbit speed button  -do not hit OK when in rabbit
- Run for a couple mins, then click play  (25 rpm) & OK
- Disconnect tubing to Autosampler & put uptake tubing into di-H₂O
- If using CCT mode, go to ‘Configurations’ & check ‘X Series Default.’ Then go to ‘tune’ and click  to reset settings

STOP HERE if someone is on after you

•Turn Off:

- Take tube out of di-H₂O and run for a couple mins
It is recommended to turn off the instrument dry.
- Turn the instrument to ‘Vacuum state’ 
- DO NOT switch the instrument to ‘Shutdown State’; this happens if you press the ‘Off’ button twice.

•Before you leave:

- ★•Release the tension on the peri-pump tubing & unhook the tubes from the tabs on the pump.
- Turn off Autosampler
- Empty waste container
- (closing Plasmalab will not shut off the ICP)
- Clean up area of any sample tubes, paper, etc.
- ★•Turn off the Argon gas or liquid tank
- ★•Log off FOM

Call me if you have any problems:

Becky: 7-7282 lab

Thank you!!